

Endomorphin-2 EIA Kit

Catalog #	EK-044-11																
Standard Size	96 wells																
Sequence	Tyr - Pro - Phe - Phe - NH ₂																
Cross Reactivity	<table><thead><tr><th><i>Peptide</i></th><th><i>%Cross-reactivity</i></th></tr></thead><tbody><tr><td>Endomorphin-2</td><td></td></tr><tr><td>Endomorphin-1</td><td></td></tr><tr><td>Orphanin FQ (Nociceptin)</td><td></td></tr><tr><td>Leu-Enkephalin</td><td></td></tr><tr><td>Dynorphin A</td><td></td></tr><tr><td>β-Endorphin (Human)</td><td></td></tr><tr><td>Met-Enkephalin</td><td></td></tr></tbody></table>	<i>Peptide</i>	<i>%Cross-reactivity</i>	Endomorphin-2		Endomorphin-1		Orphanin FQ (Nociceptin)		Leu-Enkephalin		Dynorphin A		β-Endorphin (Human)		Met-Enkephalin	
<i>Peptide</i>	<i>%Cross-reactivity</i>																
Endomorphin-2																	
Endomorphin-1																	
Orphanin FQ (Nociceptin)																	
Leu-Enkephalin																	
Dynorphin A																	
β-Endorphin (Human)																	
Met-Enkephalin																	
Range	0-100 ng/ml																
Sensitivity	Minimum Detectable Concentration = 0.1 ng/ml																
Linear Range	0.1-2.2 ng/ml																
Precision_Inter_Assay	Inter-assay error: <14%																
Precision_Intra_Assay	Intra-assay error: <5%																
Positive_Control	Rehydrate Positive Control (0.4~0.6ng/ml) with 200ul assay buffer. Add in duplicates, 50ul Standard, Controls and Samples.																

Kit Manufacturer



PHOENIX PHARMACEUTICALS, INC.

330 Beach Rd, Burlingame, CA, 94010, USA

PHONE (650) 558-8898 EMAIL info@phoenixpeptide.com

WWW.PHOENIXPEPTIDE.COM

Please see our website www.biosupply.co.uk or call us for pricing
Discounts available on bulk orders

T: 0870 760 5152

F: 0871 236 4352

E: sales@biosupply.co.uk